



Product catalogue

Contents

Protective apparel	4
Bed protection	20
Sampling & examining	24
Patient hygiene & care	30
Wound care	35
Surgical tapes	41
Adhesive tapes	42

Protective apparel

Correct use of protective apparel prevents microbes and infection from spreading from personnel and visitors to patients and vice versa. Disposable PE-apron protects better than textile protective gown, where the particles are stored in the structure of the material, and is recommended to be used when there is a risk of working clothes getting wet or soiled. PE gloves give extra protection for example when used together with washing/examination gloves.

SELEFA® vinyl, latex and nitrile gloves offer a cost efficient solution both in health care and non-medical settings

Protective apparel assortment contains shoe covers, arm covers, caps and overalls etc. Many of these products have an important role in infection control but are also very well suited to non-medical uses.

SELEFA® assortment contains both medical devices (MDD 93/42/EEC) and personal protective equipment (PPE 89/686/EEC). All medical devices in this catalogue are Class I or Class Is (sterile) and all personal protective equipment are Cat I.

The labeling follows EN 980 and EN1041 standards. All products are equipped with barcodes to help handling in the warehouse.



Comfortably in control with
SELEFA®

Comfortably in control with SELEFA® aprons

Disposable aprons made of fluid resistant polyethylene are an important part of infection control strategy. They also help protect clothing from spills and splashes during medical procedures and cleaning tasks.

All Selefa aprons are in conformity with the European Directive for Personal Protective Equipment 89/686/EEC, Cat I.

Lightweight and disposable, these aprons are also ideal for use in food preparation and in industrial and pharmaceutical applications.

Selefa polyethylene aprons are available in a variety of sizes (widths/lengths), thicknesses and colours.

The packing systems available are

- Flat packed, 210405S
- On a block hanger, 210431
- On a roll in dispenser eg 210433
- Individually folded in a dispenser eg 21043050 (Bosco)
- Individually wrapped in a dispenser eg 210424

Available in SELEFA® assortment are handy apron holders made of plexiglass.



Aprons without sleeves

Apron

REF 210415

- White, LDPE film, embossed
- Size: 72 x 118 cm
- Packed in a handy dispenser box
- Pcs/box: 100
- Pcs/carton: 500



Apron Bosco

REF 21043050

- Transparent, LDPE film, embossed
- Size: 81 x 125 cm
- Individually folded
- Packed in a handy dispenser box
- Pcs/box: 50
- Pcs/carton: 500



Apron on a roll

REF 210433

- Transparent, LDPE film, embossed
- Size: 81 x 125 cm
- Packed into a handy dispenser box
- Pcs/box: 50 on a roll
- Pcs/carton: 500



Apron on a block hanger

REF 210432

- Transparent, LDPE film, embossed
- Size: 81 x 125 cm
- Easy to dispense from hanger
- Pcs/bag: 50 on block hanger
- Pcs/carton: 500



Apron

REF 210405S

- White, LDPE film, embossed
- Size: 81 x 125 cm
- Packed into a plastic bag with cardboard support
- Pcs/bag: 100
- Pcs/carton: 500



Apron

REF 210406S

- White, LDPE film, embossed
- Size: 80 x 140 cm
- Pcs/bag: 100
- Pcs/carton: 500



Long aprons without sleeves

Apron, extra long

REF 210434

- Transparent, LDPE film, embossed
- Size: 81 x 150 cm
- Individually folded
- Packed into a handy dispenser box
- Pcs/box: 40
- Pcs/carton: 400



Apron, extra long, in a roll

REF 210435

- Transparent, LDPE film, embossed
- Size: 81 x 150 cm
- Packed into a handy dispenser box
- Pcs/box: 50
- Pcs/carton: 500



Apron, extra long

REF 210408S

- White, LDPE film, embossed
- Size: 85 x 160 cm
- Pcs/bag: 100
- Pcs/carton: 500



Apron, extra long

REF 210407S

- White, LDPE film, embossed
- Size: 90 x 160 cm
- Pcs/bag: 100
- Pcs/carton: 500



Heavy duty aprons without sleeves

Apron, extra strong

REF 211400

- White, LDPE film
- Size: 85 x 140 cm
- Pcs/bag: 100
- Pcs/carton: 500



Apron, extra strong, on a block hanger

REF 210431

- Transparent, LDPE film, embossed
- Size: 100 x 150 cm
- Extra wide
- Easy to dispense from block hanger
- Pcs/bag: 50 on block hanger
- Pcs/carton: 500



Aprons with sleeves

Apron with short sleeves

REF 210426

- Transparent, LDPE film, embossed
- Size: 98 x 120 cm
- Individually wrapped
- Packed into a handy dispenser box
- Pcs/box: 30
- Pcs/carton: 150



Apron with sleeves, thumb loop

REF 210425T

- Transparent, LDPE film, embossed
- Size: 98 x 120 cm
- Individually wrapped
- Packed into a handy dispenser box
- Pcs/box: 20
- Pcs/carton: 100



Apron with sleeves, thumb loop

REF 210425L

- Blue, LDPE film, embossed
- Size: 98 x 120 cm
- Individually wrapped
- Packed into a handy dispenser box
- Pcs/box: 20
- Pcs/carton: 100



Apron with sleeves, thumb loop

REF 210424

- Blue, LDPE film, embossed
- Individually wrapped
- With thumb loop
- Size: 94 x 140 cm
- Individually wrapped
- Pcs/unit: 1
- Pcs/carton: 100




SELEFA® PE gloves

PE gloves, Extra Strong

- Transparent, smooth
- Available in 3 sizes
- Length 34 cm
- Pcs/bag: 100
- Pcs/carton: 5000

REF 210329S
REF 210330M
REF 210331L

Small
Medium
Large




Smooth surface

PE gloves

- Transparent, smooth
- Available in 3 sizes
- Length 35 cm
- Pcs/bag: 100
- Pcs/carton: 5000

REF 210320
REF 210321
REF 210322

Small
Medium
Large




Smooth surface

PE gloves

- Transparent, embossed
- Available in 3 sizes
- Length 33 cm
- Pcs/bag: 100
- Pcs/carton: 5000

REF 21030333
REF 21030533
REF 21030433

Small
Medium
Large




Embossed surface

PE gloves

- Transparent, embossed
- Available in 4 sizes
- Length 34 cm
- Pcs/bag: 100
- Pcs/carton: 5000

REF 210298
REF 210299S
REF 210300M
REF 210301L

X-small
Small
Medium
Large




Embossed surface

SELEFA® PE gloves

PE gloves

- Transparent, embossed
- Available in 3 sizes
- Length 34 cm
- Packed in handy dispenser box
- Pcs/box: 400
- Pcs/carton: 4000

REF 210299A Small
REF 210300A Medium
REF 210301A Large




Embossed surface

PE gloves, Extra Long

- Transparent, embossed
- One-size
- Length 53 cm
- Pcs/box: 100
- Pcs/carton: 2500

REF 210312 Medium




X-long, embossed surface

PE gloves, Extra Strong

- Rose, embossed
- Available in 3 sizes
- Length 35 cm
- Pcs/bag: 100
- Pcs/carton: 5000

REF 10204 Small
REF 10206 Medium
REF 10205 Large




Embossed surface

PE gloves, Extra Strong

- Rose, embossed
- Available in 3 sizes
- Length 34 cm
- Pcs/bag: 100
- Pcs/carton: 5000

REF 210332S Small
REF 210333M Medium
REF 210334L Large



Smooth surface

Powder-free paper backed copolymer gloves

Powder-free paper backed copolymer gloves, unsterile

- Available in 3 sizes
- Length 28 cm
- Transparent
- Smooth surface
- Pcs/dispenser: 100
- Pcs/carton: 1000

REF 210319S Small
REF 210319M Medium
REF 210319L Large



Powder-free paper backed copolymer gloves, sterile

- Available in 3 sizes
- Length 28 cm
- Transparent
- Smooth surface
- Gamma-sterilized
- Pcs/unit pack: 1
- Pcs/dispenser: 100
- Pcs/carton: 2000

REF 210318S Small
REF 210318M Medium
REF 210318L Large



SELEFA® cotton gloves

SELEFA® Extra

- Natural white, strong
- Available in 3 sizes
- Length 30 cm
- Pcs/bag: 40
- Pcs/carton: 2000

REF 222884L S
REF 222885L M
REF 222886L L



SELEFA® Special

- White
- Available in 3 sizes
- Length 30 cm
- Pcs/bag: 40
- Pcs/carton: 2000

REF 222889L S
REF 222890L M
REF 222891L L



SELEFA® Standard

- White
- Available in 2 sizes
- Length 25 cm
- Pcs/bag: 40
- Pcs/carton: 2000

REF 222890 M
REF 222891 L

SELEFA® Consumer

- Natural white
- Available in 3 sizes
- Length 30 cm
- Pcs/box: 10
- Pcs/carton: 1600

REF 222895ADAS S
REF 222895ADAM M
REF 222895ADAL L



SELEFA® Basic

- Natural white
- Available in 3 sizes
- Length 30 cm
- Pcs/bag: 40
- Pcs/carton: 2000

REF 222895S S
REF 222895M M
REF 222895L L



Knitted cotton gloves with grip points

- Available in white and black
- Available in 4 sizes
- Length 22 cm
- Pairs/bag: 12
- Pairs/carton: 600
(50 bags of 12 pairs)

White:
REF 222901 S
REF 222902 M
REF 222903 L
REF 222904 XL

Black:
REF 222911 S
REF 222912 M
REF 222913 L
REF 222914 XL



SELEFA® vinyl gloves

Vinyl

Powder-free SELEFA® vinyl glove offers practical protection for many applications. The 100 % synthetic glove is somewhat softer than the standard vinyl and has a good feel and performance. It is a good alternative for those suffering from Type I allergies.

Features

- DOP-free
- No latex proteins, good alternative for those suffering from Type I allergies
- Suitable for food handling, except for fatty food
- Practical barrier for wide range of applications
- Good fit, feel and performance
- AQL 1,5

Quality standards

- Complies with EN 455



Technical information

Type	Powder-free, unsterile
Material	Vinyl
Colour	Transparent
Model	Ambidextrous, with rolled rim
Storage	Protect from heat, humidity, strong light and ozone
Size/overall length as per EN 455-2	S, M and L 240 mm
Impermeability as per EN 455-1	AQL 1,5
Durability, in original package if stored as per din 7716, ISO 2230	Min. 5 years

- Pcs/dispenser: 100
- Pcs/carton: 1000

REF 210283S S
REF 210283M M
REF 210283L L

SELEFA® latex gloves

Latex, smooth

Latex powdered examination gloves offer comfort and protection, a good combination of strength and elasticity. They perform well in heavy-duty situations and fit comfortably like a second skin. The only downside to natural rubber latex examination gloves is that being a protein based compound they may cause allergic reactions and irritation in some people – including anaphylactic reaction.

Features

- Beaded cuff for easy donning and prevention of roll back
- Low level of extractable proteins and chemical residue
- Food safe, suitable for handling all foodstuffs
- AQL 1,5

Quality standards

- Complies with EN 455



Technical information

Type	Powdered, unsterile
Material	Latex
Colour	Natural white
Model	Ambidextrous, with rolled rim
Storage	Protect from heat, humidity, strong light and ozone
Size/overall length as per EN 455-2	XS, S, M, L and XL 240 mm
Impermeability as per EN 455-1	AQL 1,5
Durability, in original package if stored as per din 7716, ISO 2230	Min. 5 years

- Pcs/dispenser: 100
- Pcs/carton: 1000

REF 210277XS	XS
REF 210277S	S
REF 210277M	M
REF 210277L	L
REF 210277XL	XL

SELEFA® nitrile gloves

Nitrile general use: blue and white

Comfortable, soft flexible powder-free nitrile gloves provide added safety in many applications. Designed with a special nitrile formulation they feel and fit like latex and allow full range of motion and excellent flexibility to minimise stress and fatigue. They do not contain natural rubber latex and are excellent alternative for those suffering from Type I allergies.

Features

- 100% latex free – made from durable, protein and powder-free nitrile, eliminates the Type 1 allergic reaction associated with natural rubber latex
- Highly elastic and soft
- Rolled rim to facilitate easy donning
- Designed to give a natural rubber-like feel
- Textured in finger tips for enhanced grip in all situations
- Food safe, suitable for handling fatty foods incl. cooking, olive and cod liver oil
- Protects against a wide variety of chemicals
- AQL 1.5

Quality standards

- Complies with EN 455:1–4 and EN 374:2

Technical information

Type	Powder-free, unsterile
Material	Nitrile butadiene rubber (NBR)
Colour	Blue/white
Model	Ambidextrous, with rolled rim
Storage	Protect from heat, humidity, strong light and ozone
Size/overall length as per EN 455-2	XS, S, M, L and XL 240 mm
Impermeability as per EN 455-1	AQL 1,5
Durability, in original package if stored as per din 7716, ISO 2230	3 years



Blue:

REF 210255XS	XS	Blue: REF 210253XS	XS
REF 210255S	S	REF 210253S	S
REF 210255M	M	REF 210253M	M
REF 210255L	L	REF 213253L	L
REF 210255XL	XL	REF 213253XL	XL

Pcs/dispenser: 100
Pcs/carton: 1000

Pcs/dispenser: 200
Pcs/carton: 2000

White:

210251XS	XS	White: 210252XS	XS
210251S	S	210252S	S
210251M	M	210252M	M
210251L	L	210252L	L
210251XL	XL	210252XL	XL

Pcs/dispenser: 100
Pcs/carton: 1000

Pcs/dispenser: 200
Pcs/carton: 2000

SELEFA® wall dispensers

Wall dispenser

Acrylic box dispenser for 1, 2 or 3 boxes of examination gloves (100 pcs/box or 150 pcs/box) or aprons.

- Pcs/unit pack: 1
- Pcs/transport carton: 1 x 10

Mini	REF 215051
Medi	REF 215052
Maxi	REF 215053

Suited for pcs of boxes

1 box

2 boxes

3 boxes



Cleaning: Mild soap solution and water or alternatively 70 % ethanol or similar.

Shoe cover

REF 211541

- Unsterile, PE film, blue
- Size: 41 x 15 cm
- Pcs/unit pack: 10
- Pcs/box: 100
- Pcs/carton: 2000



Shoe cover

REF 211545

- Unsterile, PE film, blue
- Size: 45 x 15 cm
- Pcs/unit pack: 10
- Pcs/box: 100
- Pcs/carton: 2000



Shoe cover

REF 211540

- Unsterile, PE film, white, heavy-duty
- Size: 40 x 15 cm
- Pcs/box: 50
- Pcs/carton: 1000



Arm cover, elasticated cuff

REF 210750

- Unsterile, PE film, white
- Size: 42 x 22 cm
- Pcs/unit pack: 10
- Pcs/box: 100
- Pcs/carton: 2000



Arm cover, knitted cuff

REF 210745

- Unsterile, PP nonwoven/PE film, nylon cuff, yellow
- Size: 21 x 43 cm
- Pcs/unit pack: 2
- Pcs/box: 50
- Pcs/carton: 400



SELEFA® Protective apparel

Visitor coat

REF 210800

- Unsterile, PE film, white
- Size: 125 cm
- Pcs/unit pack: 1
- Pcs/box: 12
- Pcs/carton: 600



Lab coat, velcro closure, 2 pockets, without cuff

REF 210801

- Unsterile, PP nonwoven, white
- Size: 118 cm
- Pcs/unit pack: 1
- Pcs/carton: 50



Beard cover

REF 210855

- Unsterile, PP nonwoven, white
- Size: 48 x 24 cm
- Pcs/box: 100
- Pcs/carton: 1000



Face mask 2-ply, paper

REF 210451

- Unsterile, white
- Pcs/box: 100
- Pcs/carton: 10000



SELEFA® Protective apparel

Peaked cap

- Unsterile, PP nonwoven, cardboard peak
- White/white/blue peak
- Pcs/box: 100
- Pcs/carton: 1000

REF 210763
REF 210767

M (50 cm)
L (55 cm)



Peaked cap

- Unsterile, PP nonwoven, cardboard peak
- White/white/red peak
- Pcs/box: 100
- Pcs/carton: 1000

REF 210764
REF 210768

M (50 cm)
L (55 cm)



Snood peak with hairnet

REF 210760

- Unsterile, PP nonwoven
- White
- Size: XL
- Pcs/box: 100
- Pcs/carton: 1000



Overall with hood

- Unsterile, spunbond PP nonwoven
- White
- Pcs/unit pack: 1
- Pcs/carton: 50

REF 210810
REF 210811
REF 210812

M
L
XL



Bed protection


Spreading of infection is prevented by protecting the patient mattress with a plastic mattress cover. Protection is especially important if the patient is incontinent or has exuding wounds. It is also good to protect the pillow from sweat and other secretions with a PE pillow cover. Isolation of patient bed is completed with draw sheet. OneMed assortment contains a large selection of draw sheets for different needs, also a polyethylene version.



Draw sheet Special, threadreinforced

REF 500334

- White, unsterile
- Tissue/PE laminate, fluid barriers at long edges, reinforcement threads polyester
- Size: 80 x 140 cm
- Pcs/box: 25
- Pcs/carton: 100



Draw sheet

REF 500336

- White/blue, unsterile
- Tissue/PE laminate, fluid barriers at long edges
- Size: 80 x 140 cm
- Pcs/carton: 100



Draw sheet Super, threadreinforced

REF 500337

- White, unsterile
- Tissue/tissue/PE laminate, fluid barriers at long edges, reinforcement threads polyester
- Size: 80 x 140 cm
- Pcs/carton: 100



Bed protection, roll

REF 30680

- Unsterile, PE film
- Size: 90 x 130 cm
- Pcs/box: 50
- Pcs/carton: 500



Mattress cover

REF 030660

- Blue, unsterile, CPE film
- Size: 210 x 90 x 2 cm
- Pcs/unit pack: 1
- Pcs/box: 10
- Pcs/carton: 100



Mattress cover

REF 030655

- Blue, unsterile, PE film
- Size: 100 x 240 cm
- Pcs/box: 20
- Pcs/carton: 200



Mattress cover

REF 30670

- Transparent, unsterile, PE film
- Size: 75 x 160 cm
- Pcs/carton: 100

Mattress cover

REF 30665

- Transparent, unsterile, PE film
- Size: 95 x 230 cm
- Pcs/carton: 100

Strether sheet

REF 30780

- White, unsterile, viscose/polyester/cellulose nonwoven
- Size: 70 x 200 cm
- Pcs/carton: 100



Examination table crepe

REF 500510

- Natural white, unsterile
- Size: 50 x 200 cm
- Pcs/carton: 2 rolls



Pillow case

REF 30750

- White, unsterile, viscose/polyester/cellulose nonwoven
- Size: 50 x 70 cm
- Pcs/carton: 100



Pillow case

REF 030760

- Rose, unsterile, PE film
- Size: 60 x 77 cm
- Pcs/box: 100
- Pcs/carton: 1000



SELEFA® Bed protection

Pillow case

REF 030740

- Blue, unsterile, PE film
- Size: 40 x 50 cm
- Pcs/box: 100
- Pcs/carton: 1000



Pillow bag, long side open

REF 30751

- Transparent unsterile, PE film
- Size: 50 x 70 cm
- Pcs/box: 50
- Pcs/carton: 1000



Pillow bag, short side open

REF 30752

- Light green, unsterile, PE film
- Size: 50 x 75 cm
- Pcs/carton: 100
- Pcs/carton: 1500



Bed sheet

REF 30755

- Unsterile
- White, viscose/polyester/cellulose nonwoven
- Size: 140 x 240 cm
- Pcs/carton: 50



Lifting sheet

REF 30790

- Unsterile
- White PP nonwoven + transparent PE cast film
- Size: 140 x 240 cm
- Pcs/carton: 25



Sampling & examining

The 2-layer injection swab is very absorbent. It is available both loose packed in box of 400 pcs and perforated in a roll. Each transport carton of roll-packed injection swabs contains one easy-to-use cardboard dispenser. For quick skin cleansing the swab assortment contains also an isopropyl alcohol skin cleansing swab.

The birch tree wooden spatulas – tongue depressor and cervical spatula – are packed into handy dispensers.

Two- or three-layer protective drapes are suitable for many uses. They easily provide a clean surface for sampling or for other procedures requiring clean and isolating surface.

The sampling & examining products also contains ear funnels (fitting on Heine otoscope) and a wide assortment of applicators.



Protective drape

REF 500505

- White/blue, unsterile
- Tissue/PE film laminate
- Size: 37 x 40 cm
- Pcs/box: 200
- Pcs/carton: 1200



Protective drape

REF 500506

- White/blue, unsterile
- Tissue/PE film laminate
- Size: 37 x 40 cm
- Pcs/box: 250
- Pcs/carton: 1500



Protective drape

REF 500507

- white, unsterile
- Tissue/PE film/tissue laminate
- Size: 37 x 40 cm
- Pcs/box: 200
- Pcs/carton: 1200



Swab, perforated

REF 22271

- White, unsterile
- Soft drylaid paper
- Size: 8 x 10 cm
- Pcs/carton: 500
- Pcs/carton: 6000



Injection swab, perforated

REF 22275

- White, unsterile
- Soft drylaid paper
- Size: 4 x 5 cm
- Pcs/carton: 1000
- Pcs/carton: 12000



Injection swab, loose

REF 22276

- White, unsterile
- Soft drylaid paper
- Size: 4 x 5 cm
- Pcs/box: 400
- Pcs/carton: 7200



Injection swab, isopropyl alcohol 70 %

REF 222680

- White, unsterile
- PP nonwoven, isopropyl alcohol
- Size: 3 x 3 cm
- Pcs/box: 200
- Pcs/carton: 4000



Ear funnel

REF 1201024

- Black, unsterile
- Size: 4 mm
- Pcs/box: 240
- Pcs/carton: 240



Ear funnel

REF 1201152

- Grey, unsterile
- Size: 2,5 mm
- Pcs/box: 200
- Pcs/carton: 200



Ear funnel

REF 1201022

- Black, unsterile
- Size: 2,5 mm
- Pcs/box: 240
- Pcs/carton: 240



Ear funnel

REF 1201154

- Grey, unsterile
- Size: 2,5 mm
- Pcs/box: 200
- Pcs/carton: 200



Forceps

REF 210500L

- White, unsterile
- Pcs/box: 50
- Pcs/carton: 300



SELEFA® Sampling & examining

Cotton tipped applicator, single tip

REF 120783


- Unsterile
- Size: 15 cm
- Pcs/unit pack: 100
- Pcs/box: 2400
- Pcs/carton: 24000



Cotton tipped applicator, double tip

REF 120795

- Unsterile
- Size: 8 cm
- Pcs/unit pack: 100
- Pcs/box: 2400
- Pcs/carton: 24000



Cotton tipped applicator

REF 120788


- Sterile
- Size: 15 cm
- Pcs/unit pack: 2
- Pcs/box: 200
- Pcs/carton: 2000



Cotton tipped applicator

REF 120790

- Sterile
- Size: 15 cm
- Pcs/unit pack: 10
- Pcs/box: 250
- Pcs/carton: 2500



Applicator sticks

REF 120760

- Unsterile
- Size: 15 cm
- Pcs/box: 1000
- Pcs/carton: 50000



Tongue depressor

REF 140680

- Unsterile
- Size: 15 cm
- Pcs/box: 100
- Pcs/carton: 5000



Cervical spatula

REF 160298

- Unsterile
- Size: 18 cm
- Pcs/box: 100
- Pcs/carton: 5000



Stethoscope

REF 050947

- Unsterile
- Size: 60 cm
- Pcs/unit pack: 1
- Pcs/carton: 50



Ear tips for stethoscope

REF 509451

- Unsterile
- Size: 15 x 12 mm
- Pcs/unit pack: 2
- Pcs/carton: 1000



Absorbent cotton wool in roll

REF 221102

- Unsterile
- Size: 1 kg
- Pcs/unit pack: 1
- Pcs/carton: 16 rolls



Towel

REF 21201

- White, unsterile, drylaid paper
- Size: 40 x 40 cm
- Pcs/box: 50
- Pcs/carton: 300



Patient hygiene & care

Traditional washing mitten is a basic product for patient's personal hygiene. SELEFA® range comprises of three different gloves: a drylaid paper version, drylaid paper/plastic film version and a thick needle punched nonwoven version. All are safe to use and comfortable for the patient. Of the same drylaid paper are made the patient wipes that can be used also for other cleaning and wiping. A thick drylaid paper towel is a good accessory for both washing and drying bed patients and for wiping and cleaning working tables and other surfaces.

Bathtub cover that is made of PE-plastic makes the bathing trolley new and personal for each patient at each bathing time.

Impermeable tie-on and tape fixation SELEFA® bibs protect the patient clothing well during meal times. The bibs are equipped with a pocket and are available in 2 sizes. In addition to traditional 2-layer model there is a 3-layer model. The materials of the soft bibs are crepe and PE-film.

Single packed toothbrushes with or without paste, nail and foot files and shower caps complete Patient hygiene & care assortment.



Bib

REF 500705

- Unsterile, laminate tissue/PE film
- Size: 37 x 60 cm
- Pcs/box: 100
- Pcs/carton: 900



Bib

REF 500706

- Unsterile, laminate tissue/PE film
- Size: 37 x 70 cm
- Pcs/box: 100
- Pcs/carton: 900



Bib, triplex

REF 500707

- Unsterile, laminate tissue/PE film/tissue
- Size: 37 x 70 cm
- Pcs/box: 100
- Pcs/carton: 600



Bib, tape fixation

REF 500708

- Unsterile, laminate tissue/PE film
- Size: 37 x 60 cm
- Pcs/box: 100
- Pcs/carton: 900

Washglove

REF 21101

- Unsterile, airlaid paper
- White
- Size: 16 x 22 cm
- Pcs/box: 100
- Pcs/carton: 1000



Washglove with plastic lining

REF 21102

- Unsterile, airlaid paper/PE film
- White
- Size: 16 x 22 cm
- Pcs/box: 100
- Pcs/carton: 1000



Washglove extra soft

REF 21103

- Unsterile, viscose/polyester nonwoven
- White
- Size: 16 x 22 cm
- Pcs/box: 100
- Pcs/carton: 1000



Patient wipe

REF 21104

- Unsterile, drylaid paper
- White
- Size: 20 x 20 cm
- Pcs/box: 100
- Pcs/carton: 1400



Napkin

REF 210101

- Unsterile, tissue
- Size: 33 x 33 cm
- Pcs/box: 100
- Pcs/carton: 4600



Cleaning wipe

REF 21202

- Unsterile, nonwoven, viscose/polyester/cotton/cellulose fibres
- Blue
- Size: 50 x 60 cm
- Pcs/carton: 400



Toothbrush with paste

REF 210645

- Unsterile
- Pcs/pack: 1
- Pcs/box: 100
- Pcs/carton: 1000



Toothbrush without paste

REF 210640

- Unsterile
- Pcs/pack: 1
- Pcs/box: 100
- Pcs/carton: 1000



Foot file

REF 210628

- Unsterile
- Pcs/box: 100
- Pcs/carton: 2500



Nail file

REF 210625

- Unsterile
- Pcs/box: 100
- Pcs/carton: 5000



Shower cap

REF 030635

- Unsterile, PE film
- Pcs/box: 100
- Pcs/carton: 2000

Flexible drinking straw

REF 030344

- Unsterile
- Size: 21 cm
- Pcs/unit pack: 1
- Pcs/box: 200
- Pcs/carton: 5000



Triangular sling

REF 222415

- Unsterile, white
- Size: 95 x 95 x 135 cm
- Pcs/unit pack: 1
- Pcs/carton: 50



Triangular sling

REF 222410

- Unsterile, natural white
- Size: 90 x 90 x 130 cm
- Pcs/carton: 600

Basin bag

REF 210620S

- Unsterile, transparent, PE film
- Size: 65 x 40 cm
- Pcs/box: 250
- Pcs/carton: 5000



Basin bag

REF 210621

- Unsterile, transparent, PE film
- Size: 70 x 47 cm
- Pcs/box: 50
- Pcs/carton: 1000

Bath tub cover, roll

REF 30701

- Unsterile, PE film
- Transparent
- 135 m x 135/270 cm
- Pcs/carton: 1 roll

Wound care

SELEFA® nonwoven swabs and sponges are cost-effective alternatives to traditional cotton gauze. They feature superior absorption capacity, softness and low fibre release. The raw material is high quality viscose/polyester nonwoven.

SELEFA® gauze swabs, sponges and tamponades are available in large number of sizes and packaging configurations. They are made of 17-thread gauze that complies with the stipulations of the European Pharmacopoeia.

SELEFA® products are packed into easy-to-use packaging.



REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
-----	---------------------	------	---------------	---------------	------------

Gauze swabs, unsterile, 8-ply packed in paper bag

22219513	Gauze swabs	5 x 5 cm	-	100	1000
22219524	Gauze swabs	7,5 x 7,5 cm	-	100	1000
22219535	Gauze swabs	10 x 10 cm	-	100	1000
22219546	Gauze swabs	10 x 20 cm	-	100	2000

22219513



Gauze swabs, unsterile, 12-ply packed in paper bag

22219512	Gauze swabs	5 x 5 cm	-	100	1000
22219523	Gauze swabs	7,5 x 7,5 cm	-	100	1000
22219534	Gauze swabs	10 x 10 cm	-	100	1000
22219545	Gauze swabs	10 x 20 cm	-	100	1000

Gauze swabs, presterilized, unsterile, 8-ply, packed in dispenser

22219511	Gauze swabs	5 x 5 cm	-	120	10800
22219522	Gauze swabs	7,5 x 7,5 cm	-	120	7200
22219533	Gauze swabs	10 x 10 cm	-	120	7200

22219533



Gauze swabs, presterilized, unsterile, 12-ply, packed in dispenser

22219544S	Gauze swabs	10 x 20 cm	-	50	1000
------------------	-------------	------------	---	----	------

Gauze swabs, presterilized, unsterile, 16-ply, packed in dispenser

22219555	Gauze swabs	10 x 10 cm	-	60	3600
-----------------	-------------	------------	---	----	------

222250T



280500

Triangular and flat gauze sponges, presterilized, unsterile

222250T	Triangular gauze sponges	30 x 30 x 30 mm	100	-	10000
280500	Flat gauze sponges	8 x 8 mm	100	1000	10000
222260S	Flat gauze sponges	30 x 30 mm	100	100	10000

222260S



Gauze squares, presterilized, unsterile

22217815	Gauze squares	15 x 15 cm	-	500	10000
22217820	Gauze squares	20 x 20 cm	-	500	10000
22217830	Gauze squares	30 x 30 cm	-	500	5000
22217840	Gauze squares	40 x 40 cm	-	500	4000
22217850	Gauze squares	50 x 50 cm	-	250	2500

22217820



REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
-----	---------------------	------	---------------	---------------	------------

Gauze round sponges, unsterile

222253	Gauze round sponges	Ø 30 mm	-	100	2000
222256	Gauze round sponges	Ø 70 mm	-	100	1000

Gauze round sponges, presterilized, unsterile

222250S	Gauze round sponges	Ø 30 mm	-	200	3600
222250M	Gauze round sponges	Ø 40 mm	-	150	2700
222250L	Gauze round sponges	Ø 50 mm	-	100	1800
222250XL	Gauze round sponges	Ø 65 mm	-	100	1400

222250L



2222722



2222745

Gauze swabs, sterile, 8-ply, 2 or 5 pcs/bag

2222702	Gauze swabs	5 x 5 cm	2	120	2400
2222722	Gauze swabs	7,5 x 7,5 cm	2	120	2400
2222742	Gauze swabs	10 x 10 cm	2	120	2160
2222705	Gauze swabs	5 x 5 cm	5	125	2500
2222725	Gauze swabs	7,5 x 7,5 cm	5	125	2500
2222745	Gauze swabs	10 x 10 cm	5	125	2250



Gauze swabs, sterile, 12-ply, 2 or 5 pcs/bag

222291	Gauze swabs	5 x 5 cm	2	80	1440
222290	Gauze swabs	7,5 x 7,5 cm	2	80	1600
222289	Gauze swabs	10 x 10 cm	2	80	1280
2222762	Gauze swabs	10 x 20 cm	2	40	720
222296	Gauze swabs	5 x 5 cm	5	80	1600
222295	Gauze swabs	7,5 x 7,5 cm	5	80	1600
222294	Gauze swabs	10 x 10 cm	5	80	1280
2222765	Gauze swabs	10 x 20 cm	5	40	720

222295



REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
-----	---------------------	------	---------------	---------------	------------

Gauze swabs, sterile, 8-ply, special sizes

222298	Gauze swabs	20 x 80 cm	5	25	150
222299	Gauze swabs	10 x 60 cm	5	65	390
222300	Gauze swabs	10 x 10 cm	25	150	2400

Gauze, round, triangular and flat sponges, sterile

22225030	Round gauze sponges	Ø 30 mm	5	75	1500
22225040	Round gauze sponges	Ø 40 mm	5	75	1350
22225050	Round gauze sponges	Ø 50 mm	5	50	900
222254	Round gauze sponges	Ø 70 mm	5	100	800
222250TS	Triangular gauze sponges	30 x 30 x 30 mm	5	250	5000
22226030	Flat gauze sponges	30 x 30 mm	5	250	5000



Nonwoven swabs, unsterile, 4-ply, packed in paper bag

22319512	Nonwoven swabs	5 x 5 cm	-	100	6000
22319523	Nonwoven swabs	7,5 x 7,5 cm	-	100	7500
22319534	Nonwoven swabs	10 x 10 cm	-	100	5400
22319545	Nonwoven swabs	10 x 20 cm	-	100	2700



Nonwoven swabs, presterilized, unsterile, 4-ply, packed in box

22319511	Nonwoven swabs	5 x 5 cm	-	150	13500
22319522	Nonwoven swabs	7,5 x 7,5 cm	-	150	9000
22319533	Nonwoven swabs	10 x 10 cm	-	150	9000
22319544	Nonwoven swabs	10 x 20 cm	-	75	1500



REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
Nonwoven round sponges, presterilized, unsterile					
223250S	Nonwoven round sponges	Ø 30 mm	-	200	3600
223250M	Nonwoven round sponges	Ø 40 mm	-	150	2700
223250L	Nonwoven round sponges	Ø 50 mm	-	100	1800
223250T	Triangular nonwoven sponges	30 x 30 x 30 mm	100	1000	10000

223250S



223250T



Nonwoven swabs, sterile, 4-ply, packed 2 or 5 pcs/bag

2232702	Nonwoven swabs	5 x 5 cm	2	140	2800
2232722	Nonwoven swabs	7,5 x 7,5 cm	2	140	2800
2232742	Nonwoven swabs	10 x 10 cm	2	140	2520
2232762	Nonwoven swabs	10 x 20 cm	2	70	1260
2232705	Nonwoven swabs	5 x 5 cm	5	150	3000
2232725	Nonwoven swabs	7,5 x 7,5 cm	5	150	3000
2232745	Nonwoven swabs	10 x 10 cm	5	150	2700
2232765	Nonwoven swabs	10 x 20 cm	5	40	720

2232742



2232725



Nonwoven round sponges, sterile, packed 5 pcs/bag

22325030	Round nonwoven sponges	Ø 30 mm	5	75	1500
22325040	Round nonwoven sponges	Ø 40 mm	5	75	1350
22325050	Round nonwoven sponges	Ø 50 mm	5	50	900
223259	Round nonwoven sponges	Ø 70 mm	5	50	400
223250TS	Triangular nonwoven sponges	30 x 30 x 30 mm	5	250	5000

22325040



223250TS



Nonwoven round sponges, unsterile

223253	Nonwoven round sponges	Ø 30 mm	-	100	2000
223256	Nonwoven round sponges	Ø 70 mm	-	100	1000

REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
Tamponing bandages, sterile					
222463	Tamponing bandage	1 x 50 cm	1	10	120
222464	Tamponing bandage	1 x 100 cm	1	10	120
222465	Tamponing bandage	1 x 200 cm	1	10	120
222466	Tamponing bandage	1 x 500 cm	1	10	120
222467	Tamponing bandage	2 x 50 cm	1	10	120
222468	Tamponing bandage	2 x 100 cm	1	10	120
222469	Tamponing bandage	2 x 200 cm	1	10	120
222470	Tamponing bandage	2 x 500 cm	1	10	120
222471	Tamponing bandage	5 x 50 cm	1	10	120
222472	Tamponing bandage	5 x 100 cm	1	10	120
222473	Tamponing bandage	5 x 200 cm	1	10	120
222474	Tamponing bandage	5 x 300 cm	1	10	120
222475	Tamponing bandage	5 x 500 cm	1	10	120
222476	Tamponing bandage	10 x 100 cm	1	10	120
222477	Tamponing bandage	10 x 300 cm	1	10	120



222465



222473



Tamponing bandage, x-ray detectable, sterile

226055SX	Tamponing bandage	Ø 55 mm	1	25	100
-----------------	-------------------	---------	---	----	-----

226055SX




SELEFA® Surgical tapes

Surgical tape, nonwoven with holder

A latex-free, breathable nonwoven tape.
Gentle to the skin yet adheres well.
An economical, general purpose surgical tape.

REF 2405021

- White
- Size: 1,25 cm x 9,1 m
- Pcs/dispenser: 24
- Pcs/carton: 720



REF 2405022

- White
- Size: 2,5 cm x 9,1 m
- Pcs/dispenser: 12
- Pcs/carton: 360



REF 2405023

- White
- Size: 5 cm x 9,1 m
- Pcs/dispenser: 6
- Pcs/carton: 180

Surgical tape, nonwoven

A latex-free, breathable nonwoven tape.
Gentle to the skin yet adheres well.
An economical, general purpose surgical tape.

REF 2402021

- White
- Size: 1,25 cm x 9,1 m
- Pcs/dispenser: 24
- Pcs/carton: 720

REF 2402022

- White
- Size: 2,5 cm x 9,1 m
- Pcs/dispenser: 12
- Pcs/carton: 360



REF 2402023

- White
- Size: 5 cm x 9,1 m
- Pcs/dispenser: 6
- Pcs/carton: 180

Surgical tape, Silk

Latex-free, white, silk like high strength surgical tape.
Bi-directional tearing for easy tailoring.

REF 2402011

- White
- Size: 1,25 cm x 9,1 m
- Pcs/dispenser: 24
- Pcs/carton: 720

REF 2402012

- White
- Size: 2,5 cm x 9,1 m
- Pcs/dispenser: 12
- Pcs/carton: 360

REF 2402013

- White
- Size: 5 cm x 9,1 m
- Pcs/dispenser: 6
- Pcs/carton: 180



Surgical tape, PE-film

A transparent medical pe-tape.
Perforated, that tears easily.

REF 2402031

- Transparent
- Size: 1,25 cm x 9,1 m
- Pcs/dispenser: 24
- Pcs/carton: 720

REF 2402032

- Transparent
- Size: 2,5 cm x 9,1 m
- Pcs/dispenser: 12
- Pcs/carton: 360

REF 2402033

- Transparent
- Size: 5 cm x 9,1 m
- Pcs/dispenser: 6
- Pcs/carton: 180



Injection plasters

Adhesive Plaster, Nonwoven

REF 240100F

- Skin
- Size: 24 x 70 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200

REF 240101F

- Skin
- Size: 10 x 40 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200

Adhesive Plaster, Plastic

REF 240100FT

- Transparent
- Size: 24 x 70 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200

REF 240101FT

- Transparent
- Size: 10 x 40 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200



Round plasters

Adhesive Plaster, Nonwoven

REF 240102F

- Skin
- Size: 38 x 38 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200



REF 240103F

- Skin
- Size: 20 x 40 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200



Adhesive Plaster, Plastic

REF 240102FT

- Transparent
- Size: 38 x 38 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200



REF 240103FT

- Transparent
- Size: 20 x 40 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200

Childrens plasters

Adhesive Plaster, Plastic "Pinocchio"

REF PHA20PML

- Motive
- Size: 20 x 70 mm
- Pcs/dispenser: 240

