

3M™ Littmann® Stethoscopes Product Comparison Chart

An investment you make in yourself.

Littmann stethoscopes are made for professionals who are dedicated to achieving the best patient outcomes. They offer superb acoustic performance, consistently high quality, and an outstanding value.




Lightweight		Classic				Cardiology		Electronic	
<p>6</p> <p>Acoustic Performance</p> <p>Ideal for taking blood pressure readings and making limited patient assessments.</p>		<p>7-8</p> <p>Acoustic Performance</p> <p>Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.</p>				<p>9-10</p> <p>Acoustic Performance</p> <p>Designed to meet the needs of demanding clinicians. Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation.</p>		<p>10+</p> <p>Acoustic Performance</p> <p>High-performance auscultation. Hear more of what you need to hear with amplification, ambient noise reduction, and optional connectivity.</p>	
Lightweight II S.E. Stethoscope	Select Stethoscope	Classic II Pediatric and Infant Stethoscopes	Classic II S.E. Stethoscope	Classic III™ Stethoscope	Master Classic II™ Stethoscope	Cardiology III™ Stethoscope	Master Cardiology™ Stethoscope	Electronic Stethoscope Model 3100	Electronic Stethoscope Model 3200
Double sided chestpiece 118 grams Tunable diaphragm 2 year warranty	Single sided chestpiece 120 grams Tunable diaphragm 2 year warranty	Double sided chestpiece 105, 95 grams 3 year warranty	Double sided chestpiece 135 grams Tunable diaphragm 3 year warranty	Double sided chestpiece 150 grams Two tunable diaphragms Convertible open bell 5 year warranty	Single sided chestpiece 160 grams Tunable diaphragm 3 year warranty	Double sided chestpiece 175, 180 grams Available in 22" / 56 cm and 27" / 69 cm Two tunable diaphragms Dual-lumen tubing 5 year warranty	Single sided chestpiece 175, 185 grams Available in 22" / 56 cm and 27" / 69 cm Tunable diaphragm Dual-lumen tubing 7 year warranty	185 grams 2 year warranty	Bluetooth® wireless technology Recording and playback 185 grams 2 year warranty

Tunable Technology

3M invented tunable diaphragm technology to let you hear different frequency sounds without turning over and repositioning the chestpiece.



Light contact:
Hold the tunable diaphragm side of the chestpiece with light pressure to hear low frequency sounds.



Firm pressure:
Press the chestpiece a little more firmly to hear higher frequency sounds.

Quality you can trust

All Littmann mechanical stethoscopes are designed and manufactured in the United States according to strict process and quality controls. The result: consistent quality. That means any Littmann mechanical stethoscope you pick up will provide a consistently high level of performance.

3M™ Littmann® Stethoscope Ordering Information and Tube Color Guide



		Black	Burgundy	Navy Blue	Hunter Green	Plum	Chocolate	Gray	Pine Green	Ceil Blue	Purple	Raspberry	Caribbean Blue	Orange	Seafoam Green	Red	Lilac	Royal Blue	Ocean Blue	Peach	Pearl Pink	Coral	Dark Olive Green	Sand	Lavender	Turquoise	Lemon-Lime	Emerald	Champagne	Rose Pink*	Lime Green	
Model 3200	27 in (69 cm)	3200BK27	3200BU	3200NB				3200T																								
Model 3100	27 in (69 cm)	3100BK27	3100BU	3100NB																												
Master Cardiology™	22 in (56 cm)	2159																														
	27 in (69 cm)	2160 2161 ¹ 2175 ² 2176 ⁵	2163	2164	2165 2183 ²	2167						2178											2182 ⁵									
Cardiology III™	27 in (69 cm)	3128 3131BE ¹ 3128BRS ² 3157SM ⁵ 3152RBW ⁶	3129	3130	3134	3135	3137 3137CPR ³	3136		3146		3148	3138			3140							3149	3166 ⁵	3165 ⁵	3158 ⁷	3159 ⁵	3161 ⁵				
Master Classic II™	27 in (69 cm)	2144L 2141 ¹ 2142G ⁴	2146	2147	2632				2634	2633	2143		2630																		2636 ¹	
Classic III™	27 in (69 cm)	5620 5803 ¹ 5811 ⁵	5627	5622	5624	5831	5809 ³	5621		5630		5626 5806 ⁶	5623 5807 ⁶	5629								5633		5812 ⁵		5832	5835	5839	5840		5639	5829
Classic II S.E.	28 in (71 cm)	2201 2218BE ¹	2211	2205						2813	2209	2210	2823 ⁶	2812							2822	2817										
Classic II Pediatric	28 in (71 cm)	2113		2123					2131			2122	2119 2153 ⁶	2155		2113R		2136		2154	2115											
Classic II Infant	28 in (71 cm)	2114							2132			2125 2157 ⁶	2124	2179		2114R		2156		2158	2120											
Select	28 in (71 cm)	2290	2293						2305	2301	2294	2296	2291					2298	2306	2310	2292											

Colors shown are approximate.



3M Health Care
Littmann® Brand Products
www.3m.com/healthcare
www.Littmann.com/us

3M Offices:

Albania	+35542431044	Croatia	+38512499750	Lithuania	+37052160780
Azerbaijan	+994124045050	Estonia	+3726115900	Serbia	+381112209400
Belarus	+375173727006	Kazakhstan	+77273330000	Slovakia	+421249105211
Bulgaria	+37929601911	Latvia	+37167066120	Slovenia	+38612003630


Special chestpiece / headset finishes are only available in a limited number of tubing colors.



¹ Black Edition – black finish chestpiece and eartubes.
² Brass Edition – brass finish chestpiece and eartubes.



³ Copper Edition – copper finish chestpiece and eartubes.
⁴ Gold Edition – gold finish chestpiece and eartubes.



⁶ Rainbow Edition – rainbow finish chestpiece and sunshine gold eartubes.
⁷ Rainbow Edition – rainbow finish chestpiece and brass finish eartubes.
⁵ Smoke Edition – smoke finish chestpiece and eartubes.